## Centroid computation SC4045 CONTROL FOR HIGH RESOLUTION IMAGING

João Lopes e Silva May 3, 2014 v1.0 jpedro.e.silva@gmail.com

## 1 CENTROID COMPUTATION

Devise a routine that calculates the centroids of an intensity distribution yielded by a (Shack-)Hartmann wavefront sensor. The next steps should be followed:

- 1. Given a priori knowledge of the number of lenslets/apertures that the sensor has and the ideal spot position (in the case where the phase is not aberrated), find the spots on the intensity distribution. (although this seems like a trivial task, if the aberration is very large some spots may shift a lot in the imaging plane.)
- 2. Notice that in [1, Eqs. 2.3-4], you are not told which pixel values (u, v) should you choose to perform the summation and compute the centroid. Compare different centroiding algorithms based on the analysis in [2]. Implement the one you find most suitable (after discussion with the teaching assistant).
- 3. After computing the centroid, find the expression to compute the wavefront slope based on the information in [3].
- 4. Come up with mechanisms that throw away the information of some apertures if the lighting is not enough or if a spot can't be found in the neghbourhood of its ideal position.

## REFERENCES

- [1] M. Verhaegen, "Lecture notes on control for High Resolution Imaging," May 2012.
- [2] S. Thomas, "Optimized centroid computing in a Shack-Hartmann sensor," SPIE, vol. 5490, pp. 1238 – 1246, 2004.
- [3] Spiricon, ed., Hartmann Wavefront Analyzer Tutorial. Spiricon, 2004.